

STAR 2000™



STAR Systems Enhancement Summaries

Release 18.0
October 2012

S18000121

Copyright notice

Copyright © 2012 McKesson Corporation and/or one of its subsidiaries. All Rights Reserved.

Use of this documentation and related software is governed by a license agreement. This documentation and related software contains confidential, proprietary and trade secret information of McKesson Corporation and/or one of its subsidiaries and is protected under United States and international copyright and other intellectual property laws. Use, disclosure, reproduction, modification, distribution, or storage in a retrieval system in any form or by any means is prohibited without the prior express written permission of McKesson Corporation and/or one of its subsidiaries. This documentation and related software is subject to change without notice.

Publication date

October 2012

Produced in Cork, Ireland

Product and version

STAR 2000 Release 18.0

Publication number

S18000121

Reader comments

Any comments or suggestions regarding this publication are welcomed and should be forwarded to the attention of

STAR 2000 Documentation Team
McKesson
Mail Stop ATHQ-3302
5995 Windward Parkway
Alpharetta, GA 30005

Trademarks

STAR 2000 is a trademark of McKesson Corporation and/or one of its subsidiaries.

Microsoft and Windows are registered trademarks of Microsoft Corporation in the United States and/or other countries.

All other trademarks are the property of their respective owners.

Table of Contents

About This Document **V**

 Purpose / ScopeV

 AudienceV

 Chapter OverviewV

 System RequirementsV

Chapter 1 - STAR Systems Enhancements

 OVERVIEW1-3

 Purpose1-3

About This Document

Purpose / Scope

This document contains technical and user information about new features available in Release 18.0 of STAR Systems applications. It is meant to be used in conjunction with the *STAR Systems* reference guides.

Audience

This document is intended to inform hospital personnel concerned with STAR Systems about the Release 18.0 enhancements. Such personnel may include individuals from different departments within the facility, such as the information systems, patient care or patient accounting department.

Chapter Overview

Chapter 1: STAR Systems Enhancements

This chapter provides a summary of the STIs that were completed to enhance STAR Systems functionality.

System Requirements

All Release 18.0 STAR Systems character-based enhancements require MSE Enterprise Release 18.0 or later. Any additional system requirements are provided in the enhancement summary chapter, when applicable.

Chapter 1 - STAR Systems Enhancements

OVERVIEW 1-3

 Purpose 1-3

OVERVIEW

Purpose

The following STAR Systems enhancements are part of Release 18.0:

STI	Description
M20704	<p>Functionality has been added to the HL7 interface to send an email notification when a fatal interface error occurs. The email addresses are defined for the sending application and will notify up to three email addresses or email group when a fatal error occurs. A group notification may be set up by creating a distribution list on the exchange server that looks like a real email address, such as interfacesupportgroup@hospital.com.</p> <p>The email notification contains pertinent information so that the error may be investigated and resolved. An example email is provided below. The fields are populated based on the information available in the interface audit.</p> <p>NOTE: The email notification may contain PHI (Protected Health Information) and should be intended for internal use only.</p> <p>Please contact a McKesson representative to implement this function.</p> <p>Example email:</p> <p>(ID 185) HBI - Fatal Error "New patient type cannot be blank!" Error Text: New patient type cannot be blank! Message Type: A14 Sending App: HS Message Control: 808 Field: PV1.18 Data: Name: SMITH,JANE P Acct: MRN:</p> <p>NOTE: The following STIs must be moved together: M20704, M21044, M24212, M25407, F10464</p>

■ R e a d e r C o m m e n t F o r m ■

We value your suggestions for improving our documentation. Please use this form to evaluate the *STAR Systems - Enhancement Summaries* for Release 18.0.

Topic	Poor	Fair	Good	Excellent
Organization of information	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Accuracy of information	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Completeness of information	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Clarity of information	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Amount of overview information	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Explanation of processes	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

Are there parts of this manual that could be made more helpful to you? Please explain.

Other Comments:

Thanks for your help in improving the documentation.

Your Name and Position

Hospital/Organization
Name

Telephone Number

May we contact you?

Yes or No (circle one)

Fold here

Place
Stamp
Here

STAR 2000 Documentation Team
McKesson
Mail Stop ATHQ-3302
5995 Windward Parkway
Alpharetta, GA 30005

Fold here